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手順書 (Standard Operation Procedure, SOP)

機種別資料：Navios™, Navios EX™

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International Council for Standardization in Haematology Expert Panel on Cytometry and International Society of Laboratory Hematology Task force on Platelet Counting. Platelet counting by RBC/platelet ratio method. A reference method. Am J Clin Pathol 115: 460-464, 2001 「ICSH/ISLH Plt 2001」

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- 1) International Council for Standardization in Haematology; Prepared by the expert on cytometry, Reference method for the enumeration of erythrocytes and leucocytes. Clin Lab Haemat, 1994; 16: 131-138
- 2) N. Shosaku, Extracellular vesicles and blood diseases, Int J Hemat, 2017; 105: 392-405
- 3) Herring JM et al, Microparticles in Health and Disease, J Vet Intern Med. 2013; 27: 1020-33.
- 4) S. Robet et al, Standardization of pre-analytical variables in plasma microparticle determination: results of the International Society on Thrombosis and Haemostasis SSC Collaborative workshop J Thromb Haemost, 2013; 11: 1190-1193
- 5) S. Robet et al, Standardization of platelet-derived microparticle counting using calibrated beads and a Cytomics FC500 routine flow cytometer: a first step towards multicenter studies? J Thromb Haemost, 2008; 7: 190-197.
- 6) R. Lacroix et al, Standardization of platelet-derived microparticle enumeration by flow cytometry with calibrated beads: results of the International Society on Thrombosis and Haemostasis SSC Collaborative workshop Thromb Haemost, 2010; 8: 2571-2574
- 7) ICSH and ISLH Task Force on Platelet Counting, Platelet Counting by the RBC/Platelet Ratio Method - A Reference Method, Am J Clin Pathol 2001;115:460-464
- 8) Harrison P et al, An Interlaboratory Study of a Candidate Reference Method for Platelet Counting, Am J Clin Pathol, 2001; 115: 448-459
- 9) Harrison P, Horton A, Grant D, et al Immunoplatelet counting: a proposed new reference procedure. Br J Haematol, 2000; 108: 228-235
- 10) Robert S et al, Arteriosclerosis, Thromb Vasc Biol, 2012; 32: 1054-1058
- 11) S.Cointe, et al, Standardization of microparticle enumeration across different flow cytometry platform: results of a multicenter collaborative workshop, J Thromb Haemost, 2017; 15: 187-93

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